

## CS629 Bench Top Sealer

# Qseal<sup>®</sup> • multi



---

**manufacturer**

Conroy Medical AB  
Haesthagsvaegen 14A  
SE-194 52 Upplands Vaesby  
SWEDEN

[www.conroy.se](http://www.conroy.se)

Please, contact your local  
Conroy Medical representative  
for more information.



## CS629 Qseal-multi

REF CS629:	Qseal-multi, a complete sealing system, which includes built-in Sealing head, Power cord, Rail and Operator's manual.
Type of PVC tube:	Different types and sizes of tubes up to 6.2 mm outer diameter can be sealed due to a sophisticated sensing system, which automatically adapts sealing time.
Sealing capacity:	Max. 1000 seals/hour with PVC tubes up to 5 mm outer diameter (wall thickness 0,56 mm) at 20°C (68°F). For information concerning EVA tubes, contact your local Conroy Medical representative.
Sealing time:	0.6 up to 3 sec. depending on tubing.
Intended Use:	The Qseal-multi is a fully automatic system for sealing PVC tubes connected to blood bags included either in Blood Packs or in Apheresis Disposable Sets. The sealing can be performed when the donor is still connected to the Blood Pack or the Apheresis Disposable Set.
Input Power:	100-240 V ~ 50-60 Hz
Consumption:	250 W
Power cord:	Three wire (10 A), 2.5 meter long cable with female plug, (IEC320).
RF Output:	80 W max. / 50 Ω / 40.680 MHz
Size (W x H x D)	70 x 140 x 275 mm (2.8 x 5.5 x 10.3 in)
Weight kg (lb):	1.3 kg (2.95)
Temperature:	Operating: 0 - 35°C (32 - 95°F) Storage: - 20 - 70°C (- 4 - 158°F)
Humidity:	Operating: 10 - 90% Rh (non-condensing) Storage: 10 - 90% Rh (non-condensing)
Altitude:	Operating: maximum 3000 meters (9842 feet)
In compliance with:	- EN 60601-1: 2006, General Requirements for basic safety and essential performance. - EN 60601-1-2: 2014, Collateral standards for Electromagnetic Compatibility.
Electrical safety:	Class I The Qseal-multi is used in the same environment as medical equipment (hospitals and blood banks). It must be used by highly qualified personnel.
Manufacturer according to MDD:	Conroy Medical AB Valhallavaegen 1 SE-194 63 Upplands Vaesby SWEDEN



### manufacturer

Conroy Medical AB  
Haesthagsvaegen 14A  
SE-194 52 Upplands Vaesby  
SWEDEN

### www.conroy.se

Please, contact your local  
Conroy Medical representative  
for more information.

